



ELEGANT.
INTUITIVE.
SOPHISTICATED.

StealthStation™ S8
Navigation System

Built around the way you work.

You and Stealth™. You know what you need to be efficient for your patients. For neurosurgery, spinal surgery, or ENT surgery, we have the procedural workflows that enable you to focus on what you do best—providing the best in patient care.

Medtronic

ADVANCED TECHNOLOGIES

- Two large, high definition, 27-inch touchscreen monitors provide ultimate flexibility in interacting with the system
- Two electromagnetic emitter designs, each with large tracking volume
- 1 terabyte solid-state drive—ample space to store large patient exams
- 16GB of RAM memory for fast performance and image manipulation
- Optical camera with large tracking volume provides flexibility in positioning and addresses line of sight issues

STREAMLINED WORKFLOWS

- Dual cart system design: each cart has 27-inch multi-touch monitor and independent workstation for flexibility
- Ample storage options for accessories
- Consistent user interface across all software applications

CONNECTIVITY

- Wireless connectivity to hospital and medical devices, allowing the import and export of exams from anywhere within the hospital network
- Automatic patient registration with the O-arm imaging system
- Operating system adapts to any 16:9 aspect ratio with resolutions up to full HD 1920 x 1080, for flexibility utilizing multiple screens in the OR
- Digital and analog output allows integration

CUSTOMIZED. SECURE.

- System provides both optical camera navigation and our unique electromagnetic navigation for procedural flexibility
- Each customer can define their own level of security:
 - user authentication
 - encryption
 - antivirus
 - firewall protection

Medtronic Canada

99 Hereford Street
Brampton, Ontario
L6Y 0R3
Toll free: (800) 268-5346
Telephone: (905) 460-3800

medtronic.ca

Refer to product instruction manual/package insert for instructions, warnings, precautions and contraindications. For further information, please contact Medtronic Canada at (800) 268-5346, and/or consult Medtronic's website at www.medtronic.ca.

Healthcare professionals must review the product technical manual prior to use for detailed disclosure. For information on Indications, Safety, and Warnings, call Medtronic at (800) 268-5346, or visit Medtronic's website at www.medtronic.ca.

UC201705699 EC 04.2017

© 2016 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic.

™* Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.

Medtronic